User Manual

Motion: VR Physiotherapy

Language: English
Version 1.0 - June 2025

Welcome! *Motion: VR Physiotherapy* is a virtual reality (VR) application designed to help you manage and improve lower back pain through guided, interactive exercises. With the help of a virtual trainer, you can perform targeted movements in a safe and engaging environment.

What You'll Need:

- A Meta VR headset
- VR hand controllers
- A clear play area (at least 2m x 2m)
- Comfortable clothing and a yoga mat (recommended)

Before You Begin: Make sure your play area is clear of obstacles and well-lit. Stay hydrated, and do not attempt exercises if you are feeling unwell. Always follow your physiotherapist's recommendations for duration and intensity.

What to Expect in the App: Upon launching Motion, all users are required to complete a **mandatory balance test**. This test helps assess your initial performance and movement stability. Only after completing the balance test will you be able to access the main exercise routines.

After the test, you will unlock:

- Routine selection (Stretch, Lunge, Squat etc.)
- Settings adjustment (volume, subtitles, ambient sounds)
- Difficulty settings (Easy, Medium, Hard)

Example Exercise Routine:

Routine: Mandatory Balance test

1. **NPC Demonstration** – A virtual trainer appears and performs T-pose.
2. **Your Turn** – You're prompted to copy the trainer's pose and hold it for 30 seconds.
3. **Feedback** – The system shows a text message telling you how well you did (eg. Bad, Good, Perfect)
4. **Next Pose** – You move on to a knee raise.

Each session contains 3–5 exercises and lasts around 10–15 minutes.

Safety First:

- Stop immediately if you feel pain, dizziness, or discomfort.
- Use the leave button to stop the routine.
- Keep a chair or wall nearby for extra support if needed.

For Physiotherapists: You can assign routines for your patients and monitor progress through a dedicated website. The application automatically sends session data, such as posture accuracy, time spent, and routine completion to the website. This enables you to:

- Track your patients' engagement and progress remotely
- Provide tailored advice based on real performance data
- Prescribe specific routines suited to your patient's current condition

Simply log in to the physiotherapy portal to review your patients' data and assign new sessions as needed.

Thank you for using Motion: VR Physiotherapy. We hope it helps you feel stronger, more balanced, and pain-free.

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